



1 IN 9 MEN

are diagnosed with Prostate Cancer in their lifetime.

Early screening saves lives.

Our mission at HALO Diagnostics is to offer groundbreaking diagnostic and treatment options to help more men beat prostate cancer and live fulfilling lives.

A PSA test may help detect prostate cancer earlier and has been shown to reduce death from prostate cancer by 37%.*

At HALO Diagnostics, we use multiparametric magnetic resonance imaging (mpMRI) - combined with PSA testing - for the timely and accurate detection of prostate cancer.

We are proud to partner with imaware by offering **free** access to prostate cancer screening tests.

*Rev Urol. 2010 Winter; 12(1): 66-67


PREPARING FOR A PROSTATE CANCER TEST

These are the steps to take before testing:

- ✓ There is no need to fast before taking the test.
- ✓ You need to take this test in the morning to properly capture your hormone levels.
- ✓ You need to abstain from sexual activity for 24 hours before testing.
- ✓ Return your sample to the mailbox on the same day.


HOW DOES THE TEST WORK?

The imaware™ home collection kit is easy to use, painless and requires only 3-5 drops of blood. Test requests are evaluated and if appropriate, authorized by a physician.



DRAW 3-5 DROPS OF BLOOD

Use the pain-free lancet by pressing it against the tip of your finger and press the button.



CAPTURE SAMPLE

Capture the sample into your collection device.



SECURE AND RETURN COLLECTION DEVICE

Secure the collection device, and mail it back to imaware with the provided return label.



CHECK YOUR RESULTS ONLINE

Within 5 days imaware will send an email link to view your results online and a tele-medicine professional will contact you if your PSA is elevated.

FIRST 25 REGISTRANTS GET A FREE KIT!

VISIT: www.imaware.health/halo

CODE: HALOLOSANGELES